

# implox

IMPLOX PTY. LTD.  
Medical Equipment Suppliers  
Units 22-24  
60-68 Richmond Road  
Keswick 5035  
South Australia

Telephone: (08) 8351 1455  
Facsimile: (08) 8293 7377  
email: sales@implox.com  
Web: www.implox.com



## PHILIPS HEARTSTART FRx A.E.D.

APRIL, 2020

### Philips Heartstart FRx A.E.D.

The most robust and reliable Automated External Defibrillator on the market, the Philips FRx A.E.D. has contributed to life saving successes for many years in a multitude of challenging environments.

Intuitive functionality and ease of use makes the Heartstart FRx A.E.D. the ideal companion in the early intervention of sudden cardiac arrest, especially in those potentially hostile environments.

The Philips Hearstart FRx Automated External Defibrillator is a self-maintaining device and performs daily, weekly and monthly self-checks and calibration.

Rated at **IP55**, this defibrillator is moisture resistant (drip proof), resistant to side jetting water, dust proof and crush/impact resistant.



#### *The Philips Hearstart FRx A.E.D. supplied complete with:*

##### Robust Carry Case

- 1 x Set Universal Chest Electrodes
- 1 x Disposable 4-5 year Lithium Battery
- 1 x Quick Reference Guide
- 1 x Chest Preparation Kit
- 1 x User Manual

ORDER CODE: **861305**

