

Infant scales and length measuring devices



ADE118600 Baby Scale Electronic



- Large backlight LCD display for easy reading and in dimmed surroundings
- Simple operation with four clearly marked soft touch buttons
- Automatic zero setting, tare, auto-hold, auto switch off functions
- Low battery indicator

IMPLOX PTY. LTD.
Medical Equipment Suppliers
Units 22-24
80-88 Richmond Road
Keswick 5035
South Australia

Telephone: (08) 8351 1455
Facsimile: (08) 8293 7377
email: sales@implox.com
Web: www.implox.com

- Capacity: 20kg
- Graduation: 5g
- Dimensions (Overall): 600 x 280 x 135
- Dimensions (Tray): 600x 280 x 90mm
- Weight: 3.2kg
- Power supply: Mains and battery (4 x AA) operated

ADE112600 Baby Scale Electronic Portable



- Clear 30mm display screen
- Automatic switch off & zero setting
- Tare & hold functions
- Low battery indicator
- Capacity: 20kg
- Graduation: 5g
- Dimensions (Platform): 600 x 280 x 40mm
- Dimensions (Overall): 600 x 385 x 85mm
- Weight: 2.6kg
- Power supply: mains adapter, x6 AA alkaline batteries
- Optional carry bag (pictured)

IMPLOX PTY. LTD.
Medical Equipment Suppliers
Units 22-24
80-88 Richmond Road
Keswick 5035
South Australia

Telephone: (08) 8351 1455
Facsimile: (08) 8293 7377
email: sales@implox.com
Web: www.implox.com

ADE108800 Baby & Toddler Scale Mechanical



- Robust and simple to handle
- Metal housing
- Manual zero setting & tare
- Capacity: 20kg
- Graduation: 50g < 10kg > 100g
- Dimensions: 588 x 259 343mm
- Weight: 4.5kg
- Baby tray: 550 x 265 x 80mm

ADE10027 Baby Length Measuring Mat



- Suitable for mobile use
- Fixed head piece & movable foot stop
- Skin friendly hygienic material, easy to clean
- Measuring range: 100-1000mm
- Graduation: 5mm
- Dimensions: 1260 x 280 90mm
- Weight: 516g

IMPLOX PTY. LTD.
Medical Equipment Suppliers
Units 22-24
80-88 Richmond Road
Keswick 5035
South Australia

Telephone: (08) 8351 1455
Facsimile: (08) 8293 7377
email: sales@implox.com
Web: www.implox.com

ADE10028 Baby Length Measuring Rod



- Made of anodized aluminium
- Foldable measuring tongues
- Traingle-shaped feet to support measuring rod when placed on flat surface
- Measuring range: 100-1000mm
- Graduation: 2mm
- Dimensions: 1000 x 140 x 340mm
- Weight: 700g

ADE10021 Measuring Tape (Length & Circumference) Mechanical



- Retractable band
- Measuring range: 50-1500mm
- Graduation: 1mm
- Dimensions: 55 x 80 x 20mm
- Weight: 45g

IMPLOX PTY. LTD.
Medical Equipment Suppliers
Units 22-24
80-88 Richmond Road
Keswick 5035
South Australia

Telephone: (08) 8351 1455
Facsimile: (08) 8293 7377
email: sales@implox.com
Web: www.implox.com

Column and floor scales

ADE320600/1 Column Weighing Scale Electronic



- Heavy duty electronic column scale with transport castors
- Extra wide, flat & non-slip platform
- Automatic switch off & zero setting
- Tare, hold, BMI calculation & mother/child function
- Turnable display for easy reading
- Capacity: 250kg
- Graduation: 50g < 50kg > 100g
- Dimensions (Platform): 350 x 450mm
- Dimensions (Overall): 350 x 450 x 920mm
- Weight: 7.7kg
- Power supply: mains adapter & x4 AA alkaline batteries
- Optional mechanical measuring rod (600-2100mm) available

IMPLOX PTY. LTD.
Medical Equipment Suppliers
Units 22-24
80-88 Richmond Road
Keswick 5035
South Australia

Telephone: (08) 8351 1455
Facsimile: (08) 8293 7377
email: sales@implox.com
Web: www.implox.com

ADE319660/1 Column Weighing Scale Electronic with Handrails



- Extra wide platform handrail scale with transport castors
- Automatic switch off & zero setting
- Tare, hold & BMI calculation
- User memory: 10 readings
- Capacity: 300kg / Graduation: 100g
- Dimensions (Platform): 520 x 520mm
- Dimensions (Overall): 520 x 520 x 1135mm
- Weight: 18.5kg
- Power supply: mains adapter & rechargeable (with low battery indicator) batteries

IMPLOX PTY. LTD.
Medical Equipment Suppliers
Units 22-24
80-88 Richmond Road
Keswick 5035
South Australia

Telephone: (08) 8351 1455
Facsimile: (08) 8293 7377
email: sales@implox.com
Web: www.implox.com

ADE319660/2 Column Weighing Scale Electronic with Telescopic Rod & Handrail



- Extra wide platform handrail scale with telescopic height measuring rod & transport castors
- Automatic switch off & zero setting
- Tare, hold & BMI calculation
- User memory: 10 readings
- Capacity: 300kg
- Graduation: 100g
- Measuring range: 600-2100mm
- Dimensions (Platform): 520 x 520mm
- Dimensions (Overall): 520 x 520 x 1245mm
- Weight: 20kg
- Power supply: mains adapter & rechargeable (with low battery indicator) batteries

IMPLOX PTY. LTD.
Medical Equipment Suppliers
Units 22-24
80-88 Richmond Road
Keswick 5035
South Australia

Telephone: (08) 8351 1455
Facsimile: (08) 8293 7377
email: sales@implox.com
Web: www.implox.com

ADE306800 Column Weighing Scale Dial Mechanical with Measuring Rod



- Easy to use with large clear dial and non-slip platform
- Manual zero adjustment
- Capacity: 150kg
- Graduation: 500g
- Dimensions (Platform): 270 x 370mm
- Dimensions (Overall): 280 x 665 x 1085mm
- Weight: 13.5kg
- Telescopic height measuring rod attached to column
- Measuring Range: 600-2100mm

IMPLOX PTY. LTD.
Medical Equipment Suppliers
Units 22-24
80-88 Richmond Road
Keswick 5035
South Australia

Telephone: (08) 8351 1455
Facsimile: (08) 8293 7377
email: sales@implox.com
Web: www.implox.com

ADE318800 Column Weighing Scale Beam Mechanical with Measuring Rod



- Eye-level beam scale, readable from both sides
- Supplied with telescopic height measuring rod and transport castors
- Made of steel and heavy duty iron lever system
- Manual zero adjustment
- Capacity: 200kg
- Graduation: 100g
- Measuring Range: 600-2100mm/Graduation: 1mm
- Dimensions (Platform): 275 x 375mm
- Dimensions (Overall): 275 x 560 x 1485mm
- Weight: 18kg

IMPLOX PTY. LTD.
Medical Equipment Suppliers
Units 22-24
80-88 Richmond Road
Keswick 5035
South Australia

Telephone: (08) 8351 1455
Facsimile: (08) 8293 7377
email: sales@implox.com
Web: www.implox.com

ADE320600 Floor Scale Electronic



- Strong & sturdy construction with wide platform
- Mains & battery operation
- Reversible display
- Automatic switch off & zero setting
- Tare, hold, BMI calculation & mother/child function
- Tap on start function
- Capacity: 250kg
- Graduation: 50g <50kg > 100g
- Dimensions (Platform): 380 x 305mm
- Dimensions (Overall): 400 x 350 x 44mm
- Weight: 3.8kg
- Can be calibrated

IMPLOX PTY. LTD.
Medical Equipment Suppliers
Units 22-24
80-88 Richmond Road
Keswick 5035
South Australia

Telephone: (08) 8351 1455
Facsimile: (08) 8293 7377
email: sales@implox.com
Web: www.implox.com

ADE319660 Floor Scale Bariatric Electronic



- Extra wide body platform scale
- Mains & rechargeable batteries (with low battery indicator) operation
- Automatic switch off & zero setting
- Tare, hold & BMI calculation
- User memory: 10 readings
- Capacity: 300kg
- Graduation: 100g
- Dimensions (Platform): 520 x 520mm
- Dimensions (Overall): 520 x 520 x 75mm
- Weight: 18kg
- Can be calibrated

IMPLOX PTY. LTD.
Medical Equipment Suppliers
Units 22-24
80-88 Richmond Road
Keswick 5035
South Australia

Telephone: (08) 8351 1455
Facsimile: (08) 8293 7377
email: sales@implox.com
Web: www.implox.com

ADE322600 Floor Scale Electronic



- Strong & sturdy construction with wide platform
- Mains & battery operation
- Reversible display
- Automatic switch off & zero setting
- Mother/child function
- Tap on start function
- Capacity: 250kg
- Graduation: 100g
- Dimensions (Platform): 380 x 305mm
- Dimensions (Overall): 400 x 350 x 44mm
- Weight: 3.8kg

IMPLOX PTY. LTD.
Medical Equipment Suppliers
Units 22-24
80-88 Richmond Road
Keswick 5035
South Australia

Telephone: (08) 8351 1455
Facsimile: (08) 8293 7377
email: sales@implox.com
Web: www.implox.com

ADE308800 Floor Scale Dial Mechanical



- Large non-slip platform
- Manual zero setting
- Capacity: 160kg
- Graduation: 500g
- Dimensions (Platform): 320 x 360mm
- Dimensions (Overall): 325 x 480 x 108mm
- Weight: 3.9kg

IMPLOX PTY. LTD.
Medical Equipment Suppliers
Units 22-24
80-88 Richmond Road
Keswick 5035
South Australia

Telephone: (08) 8351 1455
Facsimile: (08) 8293 7377
email: sales@implox.com
Web: www.implox.com

Measuring devices

ADE10038 Stadiometer Mechanical



- Freestanding telescopic stadiometer mounted on metal platform
- Stationary or mobile measuring
- Measuring range: 60-210cm
- Graduation: 1mm
- Dimensions: 880 x 450 x 350mm
- Weight: 6.2kg

IMPLOX PTY. LTD.
Medical Equipment Suppliers
Units 22-24
80-88 Richmond Road
Keswick 5035
South Australia

Telephone: (08) 8351 1455
Facsimile: (08) 8293 7377
email: sales@implox.com
Web: www.implox.com

ADE10042 Stadiometer Mechanical Portable



- Measuring rod is quickly taken apart into four pieces which can be stored within the platform
- Two wall distance pieces keep the rod steady and vertical to the wall without requiring fixing
- Scale can be read from both sides at the clear red indicator points
- Measuring range: 150-2100mm
- Graduation: 1mm
- Dimensions: 348 x 214 x 410mm
- Weight: 2.3kg

IMPLOX PTY. LTD.
Medical Equipment Suppliers
Units 22-24
80-88 Richmond Road
Keswick 5035
South Australia

Telephone: (08) 8351 1455
Facsimile: (08) 8293 7377
email: sales@implox.com
Web: www.implox.com

ADE10023/1 Telescopic Height Measuring Rod for Wall Mounting Mechanical



- Telescopic height measuring rod for wall mounting (& for use with electronic column scales)
- Foldable head piece for safety
- Measuring range: 60-210cm
- Graduation: 1mm
- Dimensions: 310 x 2110 x 280mm

IMPLOX PTY. LTD.
Medical Equipment Suppliers
Units 22-24
60-66 Richmond Road
Keswick 5035
South Australia

Telephone: (08) 8351 1455
Facsimile: (08) 8293 7377
email: sales@implox.com
Web: www.implox.com

ADE10017 Measuring Tape Mechanical for Wall Mounting



- Metal roll-up measuring tape for wall mounting
- Reading is taken at the indicator line in the display window
- Measuring range: 0-220cm
- Graduation: 1mm
- Dimensions: 140 x 35 x 130mm
- Weight: 180g

IMPLOX PTY. LTD.
Medical Equipment Suppliers
Units 22-24
80-88 Richmond Road
Keswick 5035
South Australia

Telephone: (08) 8351 1455
Facsimile: (08) 8293 7377
email: sales@implox.com
Web: www.implox.com

Wheelchair and chair scales

ADE500660 Platform Wheelchair Scale Electronic



- Low level platform scale with one ramp (second ramp available for purchase) & transport castors
- Powered by mains & rechargeable batteries
- Cable remote display can be wall mounted
- Low battery indicator on display
- Automatic switch off & zero setting
- Tare, hold & BMI calculation
- Capacity: 300kg
- Graduation: 100g
- Dimensions (Platform): 800 x 800mm
- Dimensions (Overall): 915 x 800 x 44mm
- Weight: 22kg

IMPLOX PTY. LTD.
Medical Equipment Suppliers
Units 22-24
80-88 Richmond Road
Keswick 5035
South Australia

Telephone: (08) 8351 1455
Facsimile: (08) 8293 7377
email: sales@implox.com
Web: www.implox.com

ADE500660/1 Platform Wheelchair Scale Electronic with Foldable Column




- Low level platform scale with one ramp (second ramp available for purchase) & transport castors
- Powered by mains & rechargeable batteries
- Foldable column housing display unit
- Low battery indicator on display
- Automatic switch off & zero setting
- Tare, hold & BMI calculation
- Capacity: 300kg
- Graduation: 100g
- Dimensions (Platform): 800 x 800mm
- Dimensions (Overall): 915 x 800 x 900mm
- Weight: 27kg

IMPLOX PTY. LTD.
Medical Equipment Suppliers
Units 22-24
80-88 Richmond Road
Keswick 5035
South Australia

Telephone: (08) 8351 1455
Facsimile: (08) 8293 7377
email: sales@implox.com
Web: www.implox.com

ADE500660/2 Platform Wheelchair Scale Electronic with Side Rails



- Low level platform scale with one ramp (second ramp available for purchase) & transport castors
- Supplied with two side rails
- Powered by mains & rechargeable batteries
- Cable remote display can be wall mounted
- Low battery indicator on display
- Automatic switch off & zero setting
- Tare, hold & BMI calculation
- Capacity: 300kg
- Graduation: 100g
- Dimensions (Platform): 800 x 800mm
- Dimensions (Overall): 915 x 800 x 200mm
- Weight: 27kg

IMPLOX PTY. LTD.
Medical Equipment Suppliers
Units 22-24
80-88 Richmond Road
Keswick 5035
South Australia

Telephone: (08) 8351 1455
Facsimile: (08) 8293 7377
email: sales@implox.com
Web: www.implox.com

ADDITIONAL SCALES AVAILABLE FOR PURCHASE

HOP2/1 Baby Scale 2 in 1



- Combination scale, removable tray top for toddler
- Direct reading on backlit display
- Battery operated (4 AA)
- Automatic switch off & zero setting
- Tare & memory functions
- Capacity: 20kg
- Graduation: 5g
- Dimensions (Tray) 254 x 571.5 x 63.5mm

IMPLOX PTY. LTD.
Medical Equipment Suppliers
Units 22-24
80-88 Richmond Road
Keswick 5035
South Australia

Telephone: (08) 8351 1455
Facsimile: (08) 8293 7377
email: sales@implox.com
Web: www.implox.com

FM650 Column Weighing Scale Electronic



Available With Printer
(FA-632PR)

- Electronic column scale with 4 transport castors
- Three options available (scale, scale + electric height measuring or scale + electric height measure + printer)
- Automatic switch off, tare, hold & BMI functions
- User memory: 25 readings
- Capacity: 250kg
- Graduation: 100g
- Dimensions (Platform): 310 x 310mm
- Weight: 8.1kg
- Power supply: mains & rechargeable batteries (with low battery indicator)
- Optional mechanical measuring rod (600-2100mm) available, graduation 1mm

IMPLOX PTY. LTD.
Medical Equipment Suppliers
Units 22-24
80-88 Richmond Road
Keswick 5035
South Australia

Telephone: (08) 8351 1455
Facsimile: (08) 8293 7377
email: sales@implox.com
Web: www.implox.com

TCS-C3 Column Weighing Scale Electronic with Height Measuring Rod (Mechanical)



- Electronic column scale with adjustable feet and integrated measuring rod 1190-2240mm, graduation 1mm
- Automatic switch off, tare, hold & BMI functions
- Capacity: 300kg
- Graduation:50g
- Dimensions (Platform): 440 x 460mm
- Power supply: mains & internal battery batteries

IMPLOX PTY. LTD.
Medical Equipment Suppliers
Units 22-24
80-88 Richmond Road
Keswick 5035
South Australia

Telephone: (08) 8351 1455
Facsimile: (08) 8293 7377
email: sales@implox.com
Web: www.implox.com

TCS-WCR Platform Wheelchair Scale Electronic with Rail



- Low level platform scale (1m x 1m) with low 20cm slope ramp both ends, transport castors and handle
- Fold up handrail
- Powered by rechargeable batteries
- 4 heavy duty load cells
- Capacity: 500kg
- Graduation: 200g
- Dimension (overall): 1.2m x 1.4m

IMPLOX PTY. LTD.
Medical Equipment Suppliers
Units 22-24
80-88 Richmond Road
Keswick 5035
South Australia

Telephone: (08) 8351 1455
Facsimile: (08) 8293 7377
email: sales@implox.com
Web: www.implox.com

B-CHAIR Chair Scale Electronic



- Wide seat chair scale with 4 braking castors
- Retractable arm rests & fold up foot rests
- Powered by rechargeable batteries
- 4 heavy duty load cells
- Capacity: 250kg
- Graduation: 50g
- Dimensions (Seat): 530 x 460mm
- Stainless steel

IMPLOX PTY. LTD.
Medical Equipment Suppliers
Units 22-24
80-88 Richmond Road
Keswick 5035
South Australia

Telephone: (08) 8351 1455
Facsimile: (08) 8293 7377
email: sales@implox.com
Web: www.implox.com